

QATCH Technologies Release Notes

	Software v2.6b40	Firmware v2.6b40
Features & Enhancements Added	Updated build info Release date 06/25/2024	Updated build info Release date 06/25/2024
	<ul style="list-style-type: none"> - Log errors to console during DEV_get_active() - 20x baseline for mode hopping - Increase sampling speeds for fast filling datas - Set 'temp_adjusted_CA_factor' to zero - ModelData throws exceptions that hard-stop - ModelData reports excessive warnings - Use a constant 72.0 for ST during Analyze - Hide user prompts and views of ST in Run Info - CA offsets based on 'concentration' removed - Incorporate new QModel for point prediction - Revise shear rate calculation from prior builds - Hide p-points outside of +/-10% of visc-region 	<ul style="list-style-type: none"> - Remove CAL blips on left of 2nd-ary devices
Bugs Fixed	- BUG FIX: Allow runs to analyze when 90 secs	- NONE
Known Issues	- NONE	- NONE