

## QATCH Technologies Release Notes

	Software v2.6b55	Firmware v2.6b53
<b>Features &amp; Enhancements Added</b>	<b>Updated build info</b> Release date 03/17/2025	<b>Updated build info</b> Release date 02/10/2025
	<ul style="list-style-type: none"> <li>- Remove fixed height on labels in Login UI</li> <li>- Update "drop effect" correction in CurveOptimizer.py</li> <li>- Merge issue #56: Adjust when audits are requested</li> <li>- Merge issue #57: Add "Run QModel Again" button</li> <li>- Merge issue #59: Adjust Drop Effect Parameters</li> <li>- Merge issue #65: Blue line does not plot like orange</li> <li>- Merge issue #70: Remove Temp CA offset support</li> <li>- Merge issue #71: Skip __MACOSX files on import</li> <li>- Hide "Temperature" column from Analyze if all NaNs</li> </ul>	- No functional changes from v2.6b46
<b>Bugs Fixed</b>	- BUG FIX: Analyze shows all NaNs from secondaries	- NONE
<b>Known Issues</b>	- Temperatures from primary needed by secondaries (Multiplex systems only, see issue <a href="#">#69</a> )	- NONE