

QATCH Technologies Release Notes

	Software v2.6r60	Firmware v2.6r60
Features & Enhancements Added	Updated build info Release date 05/19/2025	Updated build info Release date 05/16/2025
	<ul style="list-style-type: none"> - Updated QModel performance to QModel v3.15 - Fix "drop effect" corrections in Advanced Settings - Update "drop effect" in CurveOptimizer.py - Remove 'plate-config.json' parser code, for now - Remove fixed height on labels in Login UI - Hide "Temperature" column from Analyze if all NaNs - Tweak constant parameters for CurveOptimizer - Handle hiding fields that are not enabled - Resize window to fit vertically as fields show/hide - Add buffer type to run info window - Add progress bar for build extraction during updates - Merge issue #17: NEW: User Preferences support - Merge issue #23: Deleting last user account oddity - Merge issue #50: Add Caps Lock message to login - Merge issue #51: Fix IOB Error in Analyze - Merge issue #56: Adjust when audits are requested - Merge issue #57: Add "Run QModel Again" button - Merge issue #58: NEW: Real-time run fill forcaster - Merge issue #59: Adjust Drop Effect Parameters - Merge issue #61: Analyze error with empty slice - Merge issue #64: Analyze with excluded POI is weird - Merge issue #65: Blue line does not plot like orange - Merge issue #69: (Multi) Move temps to secondaries - Merge issue #70: Remove Temp CA offset support - Merge issue #71: Skip __MACOSX files on import - Merge issue #72: Apply drop-effect to High-Shear - Merge issue #74: BUG: Disable Advanced button - Merge issue #78: Move QModel button to toolbar - Merge issue #81: QModel Upgrade to Version 3 - Merge issue #84: VisQAI mode coming soon - Merge issue #86: NEW: Advanced Model Selection - Merge issue #88: Revamp Run Info fields - Merge issue #90: Adjust density with "sucrose" 	<ul style="list-style-type: none"> - Merge issue #94: Fan stays on after cooldown ends - Align fan off event with displayed countdown timer
Bugs Fixed	- NONE	- NONE
Known Issues	- NONE	- NONE