

QATCH Technologies Release Notes

	Software v2.6b61	Firmware v2.6b60
Features & Enhancements Added	Updated build info Release date 07/02/2025	Updated build info Release date 05/16/2025
	<ul style="list-style-type: none"> - Disable partial fill predictor with QModel v3 - Improvements to QModel point prediction - Improvements to drop effect correction - Add partial fill predictor to Analyze mode - Re-enable polynomial correction with new coeffs - Add drop effect PDF to analyze results ZIP - Initial support for VisQ.AI preview for demo - Add "Advanced Information" tab to Run Info 	<ul style="list-style-type: none"> - Merge issue #94: Fan stays on after cooldown ends - Align fan off event with displayed countdown timer
Bugs Fixed	- NONE	- NONE
Known Issues	- NONE	- NONE