

QATCH Technologies Release Notes

	Software v2.6b37	Firmware v2.6b37
Features & Enhancements Added	Updated build info Release date 05/27/2024	Updated build info Release date 05/27/2024
	<ul style="list-style-type: none"> - Reduce CA offsets for protein formulations to from 7 and 5 to 5 and 3, respectively - Remove "Press Stop" from Plots UI until improved - Update UUIP in targets database if it has changed - Show expected checksum while calculating actual - Remove inaccurate text from EXE launcher script - Offer to update local Python environment after updating requirements file 	<ul style="list-style-type: none"> - Disable Ethernet chip in NV features by default Use command 'EEPROM 4 1' to enable Ethernet - Initialize 'ambient' variable to NAN in global scope
Bugs Fixed	<ul style="list-style-type: none"> - BUG: Resolve requirements while honoring any constraints put on other module dependencies - BUG: SW updates not working in SW v2.6r36_exe due to improper version of 'requests' module for 'dropbox' API to work (constraint not honored) - BUG: Only call to '_dbx_connection' if it exists - BUG: Set a default value to 'res_download' if none - BUG: Only launch EXE after an update if the checksum verification passes 	- NONE
Known Issues	- NONE	- NONE