

QATCH Technologies Release Notes

| | Software v2.6b41 | Firmware v2.6b40 |
|--|--|---|
| Features & Enhancements Added | Updated build info Release date 06/28/2024 | Updated build info Release date 06/25/2024 |
| | | |
| | <ul style="list-style-type: none"> - Add MODEL parameters to constants.py - Remove unused import from ModelData - Baseline from 0.5 - 2.5 secs in ModelData - Allow models to be disabled from constants.py - Disable the new QModel predictor until it is refined - Add option to bundle XGBoost into compiled EXE - Revert existing modules to their v2.6b39 versions | <ul style="list-style-type: none"> - Remove CAL blips on left of 2nd-ary devices |
| Bugs Fixed | - NONE | - NONE |
| | | |
| Known Issues | - NONE | - NONE |
| | | |