Release Notes

QATCH Technologies Release Notes

	Software v2.6b48	Firmware v2.6b46
Features & Enhancements Added	- Improved Al support for long runs > 90 secs - Prevent out-of-order POIs 4 and 5 predictions - Improve POI 6 predictions with more adjustment - Add 3 parameters to constants for drop effect - Add drop effect option to Advanced Settings - Reload active run when turning on/off drop effect - Apply drop effect diff_factor and vector offset	Updated build info Release date 08/18/2024 - Tread PIDs A-D as secondary devices
Bugs Fixed	- NONE	- BUG FIX: Calibration sweeps were slowed in r45
Known Issues	- NONE	- NONE