

QATCH Technologies Release Notes

	Software v2.6b44	Firmware v2.6b40
Features & Enhancements Added	Updated build info Release date 07/30/2024	Updated build info Release date 06/25/2024
	<ul style="list-style-type: none"> - Use ST = 60 when concentration > 20mg/ml - Simplify limits on initial fill trendline check - Check the point_factor_limit is valid - Increase debug logging for initial fill trendline check - Add region limit to QModel predict of CH1 and CH2 - Use QModel for all points except for the first point 	<ul style="list-style-type: none"> - Remove CAL blips on left of 2nd-ary devices
Bugs Fixed	- NONE	- NONE
Known Issues	- NONE	- NONE