## Release Notes

## QATCH Technologies Release Notes

	Software v2.6b38	Firmware v2.6b37
Features & Enhancements Added	- Initial support for GitHub update engine - Add progress bar when updating resource files - Make Amplitude plot for MULTI not clip at edges - Stop MULTI measurement in FW on Stop in SW - Changes to Plots UI status label to work on MULTI - Do not adjust contact angle if 'concentration' < 1	Updated build info Release date 05/27/2024  - Disable Ethernet chip in NV features by default Use command 'EEPROM 4 1' to enable Ethernet - Initialize 'ambient' variable to NAN in global scope
Bugs Fixed	- BUG: Apply correct limits to plot2 for ch1 and 4	- NONE
Known Issues	- NONE	- NONE